

## Product Datasheet

### Mouse Glutathione peroxidase 3 ELISA Kit (orb1658392)

|                            |                                                                                                                                                                                                 |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse Glutathione peroxidase 3<br>ELISA Kit                                                                                                                                                     |
| <b>Reactivity</b>          | Mouse                                                                                                                                                                                           |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                                           |
| <b>Note</b>                | For research use only                                                                                                                                                                           |
| <b>Application notes</b>   | Detect Range: 0.156-10 ng/mL. Assay<br>Time: 3-5 h. Detection wavelength:<br>450 nm. Sample Volume: 100 T.<br>Sample Type: Serum, plasma, tissue<br>homogenates and other biological<br>fluids. |
| <b>Expiration Date</b>     | Please enquire.                                                                                                                                                                                 |