

Product Datasheet

CLEC2D antibody (orb165875)

Description	Rabbit polyclonal antibody to CLEC2D
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	CLEC2D
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	21849
Uniprot ID	Q9UHP7
NCBI	NP_037401.1 , NP_001184248.1 , NP_001184246.1 , NP_001184247.1 , NP_001004419.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

