

Product Datasheet

Human E3 ubiquitin-protein ligase HACE1 ELISA Kit (orb1658793)

Description	Human E3 ubiquitin-protein ligase HACE1 ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: less than 0.054 ng/mL. Assay Time: 3-5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: Serum,plasma,tissuehomogenatesandotherbiologicalfluids.
Expiration Date	Please enquire.