

Product Datasheet

PABPC1L2B antibody (orb165884)

Description Rabbit polyclonal antibody to PABPC1L2B (C-term).

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications WB

Immunogen Synthetic Peptide

Target PABPC1L2A

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal

MW 22799

Uniprot ID [Q5JQF8](#)

NCBI [NP_001012995.1](#), [NP_001035971.1](#)

Dilution Range WB: 1:1000

Expiration Date 12 months from date of receipt.

