

Product Datasheet

THBS3 antibody (orb165894)

Description	Rabbit polyclonal antibody to THBS3 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	THBS3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	104201
Uniprot ID	P49746
NCBI	NP_001239537.1 , NP_009043.1 , NP_001239536.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

