

Product Datasheet

POLR1A antibody (orb165897)

Description	Rabbit polyclonal antibody to POLR1A (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	POLR1A
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	194811
Uniprot ID	O95602
NCBI	NP_056240.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

