

Product Datasheet

CRELD2 antibody (orb165902)

Description	Rabbit polyclonal antibody to CRELD2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	CRELD2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	38192
Uniprot ID	Q6UXH1
NCBI	NP_001128573.1 , NP_077300.3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

