

Product Datasheet

TRPV4 antibody (orb165910)

Description	Rabbit polyclonal antibody to TRPV4 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	TRPV4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	98281
Uniprot ID	Q9HBA0
NCBI	NP_001170904.1 , NP_001170899.1 , NP_671737.1 , NP_001170902.1 , NP_067638.3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

