

Product Datasheet

Mouse CD59A glycoprotein ELISA Kit (orb1659107)

Description	Mouse CD59A glycoprotein ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 78-5000 pg/mL. Sensitivity: less than 20 pg/mL. Assay Time: 3-5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.