

Product Datasheet

TFEB antibody (orb165914)

Description	Rabbit polyclonal antibody to TFEB
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	TFEB
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	52865
Uniprot ID	P19484
NCBI	NP_001258873.1 , NP_001161299.2 , NP_001258874.1 , NP_009093.1 , NP_001258872.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

