

Product Datasheet

CATSPER2 antibody (orb165918)

Description	Rabbit polyclonal antibody to CATSPER2 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	CATSPER2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	62041
Uniprot ID	Q96P56
NCBI	NP_742093.1 , NP_473361.1 , NP_742095.1 , NP_001269238.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

