



Product Datasheet

ZNF500 antibody (orb165934)

Description	Rabbit polyclonal antibody to ZNF500 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ZNF500
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	53674
Uniprot ID	O60304
NCBI	NP_067678.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



WB analysis of K562 cell line using ZNF5...