

## Product Datasheet

### Human Protein FAM214A ELISA Kit (orb1659810)

|                            |                                                                                                                                                                                                                       |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Human Protein FAM214A ELISA Kit                                                                                                                                                                                       |
| <b>Reactivity</b>          | Human                                                                                                                                                                                                                 |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                                                                 |
| <b>Note</b>                | For research use only                                                                                                                                                                                                 |
| <b>Application notes</b>   | Detect Range: 0.78-50 ng/mL.<br>Sensitivity: less than 0.39 ng/mL. Assay Time: 3-5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids. |
| <b>Expiration Date</b>     | Please enquire.                                                                                                                                                                                                       |