

Product Datasheet

Mouse Alpha-2A adrenergic receptor ELISA Kit (orb1659859)

Description	Mouse Alpha-2A adrenergic receptor ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Assay Time: 3-5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.