

Product Datasheet

TRIX3 antibody (orb165992)

Description	Rabbit polyclonal antibody to TRIX3 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	This TRIX3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 41-67 amino acids from the N-terminal region of human TRIX3.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Uniprot ID	A6NDQ2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

