

Product Datasheet

CD200R1L antibody (orb166049)

Description	Rabbit polyclonal antibody to CD200R1L (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	CD200R1L
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	29920
Uniprot ID	Q6Q8B3
NCBI	NP_001186144.1 , NP_001008784.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

