

Product Datasheet

Human Aldehyde dehydrogenase family 16 member A1 ELISA Kit (orb1660492)

Description Human Aldehyde dehydrogenase family 16 member A1 ELISA Kit

Reactivity Human

Tested ELISA

Applications

Note For research use only

Application notes Detect Range: 0.78-50 ng/mL.
Sensitivity: less than 0.31 ng/mL.
Assay Time: 3-5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Expiration Date Please enquire.