

Product Datasheet

Human Interferon alpha-inducible protein 6 ELISA Kit (orb1660519)

Description	Human Interferon alpha-inducible protein 6 ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.2-2000 pg/mL. Sensitivity: less than 7.9 pg/mL. Assay Time: 3-5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.