

Product Datasheet

PANX2 antibody (orb166060)

Description	Rabbit polyclonal antibody to PANX2 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PANX2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	74447
Uniprot ID	Q96RD6
NCBI	NP_443071.2 , NP_001153772.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

