

Product Datasheet

CDKN2A antibody (orb166092)

Description	Rabbit polyclonal antibody to CDKN2A
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	Cdkn2a
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	19238
Uniprot ID	Q64364
NCBI	NP_034007.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

