

Product Datasheet

Cbx5 antibody (orb166093)

Description	Rabbit polyclonal antibody to Cbx5 (C-term).
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	Cbx5
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	22186
Uniprot ID	Q61686
NCBI	NP_001103686.1 , NP_001070257.1 , NP_031652.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

