

Product Datasheet

EPB41L5 antibody (orb166099)

Description	Rabbit polyclonal antibody to EPB41L5
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	EPB41L5
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	81856
Uniprot ID	Q9HCM4
NCBI	NP_001171868.1 , NP_001171867.1 , NP_001171866.1 , NP_065960.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

