

Product Datasheet

DEFA3 antibody (orb166118)

Description	Rabbit polyclonal antibody to DEFA3 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	DEFA3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	10245
Uniprot ID	P59666
NCBI	NP_005208.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

