

## Product Datasheet

### Sheep Interleukin 1 Beta (IL1b) ELISA Kit (orb1661210)

|                            |                                                                                                                                                                                                                                                         |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Sheep Interleukin 1 Beta (IL1b) ELISA Kit                                                                                                                                                                                                               |
| <b>Reactivity</b>          | Sheep                                                                                                                                                                                                                                                   |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                                                                                                   |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                                   |
| <b>Application notes</b>   | Detect Range: 31.2-2000 pg/mL.<br>Sensitivity: 9.8 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm.<br>Sample Volume: 50 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| <b>Expiration Date</b>     | Please enquire.                                                                                                                                                                                                                                         |