

Product Datasheet

Mouse Endothelin 2 (EDN2) ELISA Kit (orb1661274)

Description	Mouse Endothelin 2 (EDN2) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.1-16 ng/mL. Assay Time: 3-5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: Serum, plasma, Urine,tissue homogenates and other biological fluids.
Expiration Date	Please enquire.