

Product Datasheet

Human CHRNA7/Neuronal acetylcholine receptor subunit alpha-7 ELISA Kit (orb1661312)

Description	Human CHRNA7/Neuronal acetylcholine receptor subunit alpha-7 ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.083 ng/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: CHRNA7. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.