

Product Datasheet

Mouse B4galt1/Beta-1,4-galactosyltransferase 1 ELISA Kit (orb1661333)

Description	Mouse B4galt1/Beta-1,4-galactosyltransferase 1 ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.049 ng/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: B4galt1. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.