

Product Datasheet

Human ATF2/Cyclic AMP-dependent transcription factor ATF-2 ELISA Kit (orb1661379)

Description	Human ATF2/Cyclic AMP-dependent transcription factor ATF-2 ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 1.56-100 ng/mL. Sensitivity: 0.78 ng/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: ATF2. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.