

Product Datasheet

PMF1 antibody (orb166140)

Description	Rabbit polyclonal antibody to PMF1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PMF1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	23339
Uniprot ID	Q6P1K2
NCBI	NP_001186590.1 , NP_001186582.1 , NP_001186583.1 , NP_009152.2 , NP_001186592.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

