

## Product Datasheet

# Mouse Capn1/Calpain-1 catalytic subunit ELISA Kit (orb1661461)

|                            |                                                                                                                                                                                                                                                                   |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse Capn1/Calpain-1 catalytic subunit ELISA Kit                                                                                                                                                                                                                 |
| <b>Reactivity</b>          | Mouse                                                                                                                                                                                                                                                             |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                                                                                                             |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                                             |
| <b>Application notes</b>   | Detect Range: 0.312-20 ng/mL.<br>Sensitivity: 0.168 ng/mL. Assay Type: Sandwich. Assay Time: 3-5 h.<br>Detection wavelength: 450 nm.<br>Gene_full_name: Capn1. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids. |
| <b>Expiration Date</b>     | Please enquire.                                                                                                                                                                                                                                                   |