

Product Datasheet

Mouse Apod/Apolipoprotein D ELISA Kit (orb1661477)

Description Mouse Apod/Apolipoprotein D ELISA Kit

Reactivity Mouse

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 0.312-20 ng/mL.
Sensitivity: 0.092 ng/mL. Assay Type:
Sandwich. Assay Time: 3-5 h. Detection
wavelength: 450 nm. Gene_full_name:
Apod. Sample Volume: 100 µL. Sample
Type: Serum, plasma, tissue
homogenates and other biological
fluids.

Expiration Date Please enquire.