

Product Datasheet

Rat Cdk2/Cyclin-dependent kinase 2 ELISA Kit (orb1661561)

Description Rat Cdk2/Cyclin-dependent kinase 2
ELISA Kit

Reactivity Rat

Tested ELISA

Applications

Note For research use only

Application notes Detect Range: 78-5000 pg/mL.
Sensitivity: 39 pg/mL. Assay Type:
Sandwich. Assay Time: 3-5 h.
Detection wavelength: 450 nm.
Gene_full_name: Cdk2. Sample
Volume: 100 µL. Sample Type:
Serum, plasma, tissue homogenates
and other biological fluids.

Expiration Date Please enquire.