

Product Datasheet

Mouse C2/Complement C2 ELISA Kit (orb1661600)

Description	Mouse C2/Complement C2 ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.78-50 ng/mL. Sensitivity: 0.21 ng/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: C2. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.