

Product Datasheet

Mouse Csf1r/Macrophage colony-stimulating factor 1 receptor ELISA Kit (orb1661644)

Description	Mouse Csf1r/Macrophage colony-stimulating factor 1 receptor ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 78-5000 pg/mL. Sensitivity: 38 pg/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: Csf1r. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.