

Product Datasheet

Human CD8B/T-cell surface glycoprotein CD8 beta chain ELISA Kit (orb1661691)

Description	Human CD8B/T-cell surface glycoprotein CD8 beta chain ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 78-5000 pg/mL. Sensitivity: 24 pg/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: CD8B. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.