

Product Datasheet

PPP2R5D antibody (orb166171)

Description	Rabbit polyclonal antibody to PPP2R5D (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PPP2R5D
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	69992
Uniprot ID	Q14738
NCBI	NP_851307.1 , NP_001257405.1 , NP_851308.1 , NP_006236.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

