

Product Datasheet

Human CLEC7A/C-type lectin domain family 7 member A ELISA Kit (orb1661740)

Description	Human CLEC7A/C-type lectin domain family 7 member A ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.2-2000 pg/mL. Sensitivity: 15.6 pg/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: CLEC7A. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.