

Product Datasheet

Human CGREF1/Cell growth regulator with EF hand domain protein 1 ELISA Kit (orb1661801)

Description	Human CGREF1/Cell growth regulator with EF hand domain protein 1 ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.194 ng/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: CGREF1. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.