

Product Datasheet

Human ALOX15B/Arachidonate 15-lipoxygenase B ELISA Kit (orb1661822)

Description	Human ALOX15B/Arachidonate 15-lipoxygenase B ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 6.25-400 nmol/L. Sensitivity: 2.9 nmol/L. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: ALOX15B. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.