

Product Datasheet

Human ATF1/Cyclic AMP-dependent transcription factor ATF-1 ELISA Kit (orb1661837)

Description	Human ATF1/Cyclic AMP-dependent transcription factor ATF-1 ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 6.25-400 pg/mL. Sensitivity: 3.125 pg/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: ATF1. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.