

Product Datasheet

TBC1D3F antibody (orb166187)

Description	Rabbit polyclonal antibody to TBC1D3F
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	TBC1D3F
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	62185
Uniprot ID	A6NERO
NCBI	NP_115634.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

