

Product Datasheet

APG4B (G254) antibody (orb166190)

Description	Rabbit polyclonal antibody to APG4B (G254)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ATG4B
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	44294
Uniprot ID	Q9Y4P1
NCBI	NP_847896.1 , NP_037457.3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

