

Product Datasheet

Human ADAD1/Adenosine deaminase domain-containing protein 1 ELISA Kit (orb1661973)

Description	Human ADAD1/Adenosine deaminase domain-containing protein 1 ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 6.25-400 pg/mL. Sensitivity: 2.1 pg/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: ADAD1. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.