

Product Datasheet

Mouse Cdk5/Cyclin-dependent-like kinase 5 ELISA Kit (orb1662072)

Description	Mouse Cdk5/Cyclin-dependent-like kinase 5 ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 12.5-800 ng/mL. Sensitivity: 6.34 ng/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: Cdk5. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.