

Product Datasheet

Rat Ccnd3/G1/S-specific cyclin-D3 ELISA Kit (orb1662115)

Description	Rat Ccnd3/G1/S-specific cyclin-D3 ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 3.12-200 ng/mL. Sensitivity: 1.63 ng/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: Ccnd3. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.