

Product Datasheet

Bovine APOH/Beta-2-glycoprotein 1 ELISA Kit (orb1662177)

Description	Bovine APOH/Beta-2-glycoprotein 1 ELISA Kit
Reactivity	Bovine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 3.12-200 ng/mL. Sensitivity: 1.51 ng/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: APOH. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.