

Product Datasheet

Mouse Acat1/Acetyl-CoA acetyltransferase, mitochondrial ELISA Kit (orb1662201)

Description	Mouse Acat1/Acetyl-CoA acetyltransferase, mitochondrial ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 7.8-500 pg/mL. Sensitivity: 4.18 pg/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: Acat1. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.