

Product Datasheet

Human CDO1/Cysteine dioxygenase type 1 ELISA Kit (orb1662209)

Description	Human CDO1/Cysteine dioxygenase type 1 ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.5-32 ng/mL. Sensitivity: 0.23 ng/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: CDO1. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.