

Product Datasheet

Chicken BMP4/Bone morphogenetic protein 4 ELISA Kit (orb1662279)

Description	Chicken BMP4/Bone morphogenetic protein 4 ELISA Kit
Reactivity	Gallus
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 15.6-1000 pg/mL. Sensitivity: 5.4 pg/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: BMP4. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.