

## Product Datasheet

### Rat Protein CYR61 ELISA Kit (orb1662369)

|                            |                                                                                                                                                                                            |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rat Protein CYR61 ELISA Kit                                                                                                                                                                |
| <b>Reactivity</b>          | Rat                                                                                                                                                                                        |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                                      |
| <b>Note</b>                | For research use only                                                                                                                                                                      |
| <b>Application notes</b>   | Detect Range: 0.31 ng/mL-20.0 ng/mL.<br>Assay Time: 3-5 h. Detection wavelength: 450 nm. Sample Volume: 100 T. Sample Type: Serum, plasma, tissue homogenates and other biological fluids. |
| <b>Expiration Date</b>     | Please enquire.                                                                                                                                                                            |